


# Aspicath Thrombus Aspiration Catheter



*Biometrix Thrombus Aspiration Catheters  
provide the best solution for thrombus removal*

# Aspicath Thrombus Aspiration Catheter

*The Biometrix Thrombus Aspiration Catheter incorporates a smooth tip design to increase pushability and a unique hydrophilic coating to escalate trackability. The Aspicath's excellent pushability is also due to the innovative shaft design and the kink resistance, specifically designed for aspiration catheters.*




Aspiration Catheter Set

Code	Description
BG-A6F-ASPICATH	6F Aspiration Catheter
BG-A7F-ASPICATH	7F Aspiration Catheter

Qty/Box: 1

Technical Data	Aspicath 6F	Aspicath 7F
Range	6F	7F
Wire compatibility	0.014" (0.36mm)	0.014" (0.36mm)
Working length	145cm	145cm
Distal part flexibility/trackability	Pebax flexibility	Pebax flexibility
Proximal part push/kinking resistance	Reinforced PEEK for support	Reinforced PEEK for support
Outer diameter (O.D.)	0.067" (1.70mm)	0.078" (1.98mm)
Extraction lumen area (Proximal)	0.95mm <sup>2</sup>	1.54mm <sup>2</sup>
Extraction lumen area (Distal)	0.93mm <sup>2</sup>	1.43mm <sup>2</sup>
Aspiration speed	1.6cc/s.	2.8cc/s.



Arthesys  
RUE RENE RAZEL 4  
COSMOS BUILDING  
SACLAY 91400  
France

Tel +33 (0)1 41 11 87 77  
Fax +33 (0)1 41 11 87 70  
Email: [info@arthesys.com](mailto:info@arthesys.com)  
[www.arthesys.com](http://www.arthesys.com)



Angio\_cathMKTaspicathv.04en/18

**BIOMETRIX**  
CRITICAL CARE SOLUTIONS