



BIOMETRIX™ PRODUCT CATALOGUE

ANGIOGRAPHY & CATH-LAB


DRAINAGE & NEPHROSTOMY CATHETERS




Interventional Cardiology and Radiology healthcare professionals rely on Biometrix™ products for a complete line of Diagnostic, Therapeutic Devices and Accessories.

Through constructive dialogues with partners and customers, we are dedicated to developing effective solutions and contributing to global medical advancement.

Drainage Catheter Insertion Set




Nitinol Guide Wire with Floppy Tip




Code	Description	Qty/Box
HD-N6410	Drainage Catheter Insertion Set	5

Set Contents:


- 21Gx15cm Access Needle Two Part
- 21 Ga Stiffening Cannula
- 4Fr x 20cm Dilator
- 0.018"x 65cm Nitinol Guidewire
- 6Fx18cm Sheath with Hydrophilic Coating



Nitinol Guide Wire with Floppy Tip



Amplatz Guide Wire




Code	Description	Qty/Box
HD-N6411	Drainage Catheter Insertion Set with 0.035"x150cm Amplatz Guidewire	5
HD-N6412	Drainage Catheter Insertion Set with 0.035"x80/150cm Amplatz Guidewire	5

Set Contents:

- 21Gx15cm Access Needle Two Part
- 21 Ga Stiffening Cannula
- 0.018"x 65cm Nitinol Guidewire
- 6Fx18cm Sheath with Hydrophilic Coating
- 4Fr x 20cm Dilator
- 0.035"x 80/150cm Amplatz Guidewire

Drainage Catheter Connectors

Code	Length	Stopcock	Qty/Box
DQ-0305	5cm	-	50
DQ-0325	25cm	-	50
DQ-0430	40cm	+	100





QMD™ provides a full range of products and services to help OEMs bring medical devices to market. Specializing in silicone and other elastomerics, we develop and manufacture a broad array of high-performance single-use devices as well as precision components and subassemblies. Our end-to-end services include engineering & development, regulatory affairs, quality assurance, and supply chain support. Our global team partners with OEMs to provide tailored solutions to meet your needs.

Biometrix™

Biometrix™ offers comprehensive selection of Biometrix™ branded medical devices for endosurgery applications. Biometrix™ expertly engineers solutions and provides supply chain management for an array of clinical areas such as emergency units, operating rooms and Cath-labs.


DSP MED S.L.
QMD's Middle East Exclusive Distributor

CONTACT US >>>

DSP^{S.L.}med
Biometrix™ Products
Performance for Life

 | WWW.DSPMed.COM

 | info@dspmed.com

 | Plaza de la Solidaridad nº12, 5ª
planta 29006 Málaga
Spain