



BIOMETRIX™ PRODUCT CATALOGUE

ANGIOGRAPHY & CATH-LAB

DRAINAGE & NEPHROSTOMY CATHETERS



Interventional Cardiology and Radiology healthcare professionals rely on Biometrix™ products for a complete line of Diagnostic, Therapeutic Devices and Accessories.

Through constructive dialogues with partners and customers, we are dedicated to developing effective solutions and contributing to global medical advancement.

Guide Wires

Angiography / Diagnostics

Biometrix PTFE coated Guide Wires provide smooth consistent coverage that resist flaking and can be used for a variety of angiography procedures.



Code	Sizes	Tip Shape	Flexibility Area	Features
HR-3831	0.035"x150cm	J 3mm	Double Distal	Fixed Core, Finger Straightenable
HR-5521	0.035"x150cm	J 3mm	Single Distal	Movable Core
HR-3822	0.035"x150cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-3823	0.035"x150cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable, Firm Body
HR-3812	0.035"x150cm	Straight	Single Distal	Fixed Core
HR-3031	0.035"x180cm	J 3mm	Double Distal	Fixed Core, Finger Straightenable
HR-3922	0.035"x260cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-3912	0.035"x260cm	Straight	Single Distal	Fixed Core, Finger Straightenable
HR-4822	0.038"x150cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-4862	0.038"x150cm	Straight	Single Distal	Fixed Core
HR-4922	0.038"x260cm	J 3mm	Single Distal	Fixed Core, Finger Straightenable
HR-4912	0.038"x260cm	Straight	Single Distal	Fixed Core

Qty/Box: 50

Hydrophilic

Hydrophilic Guide wires are designed to provide superior lubricity for enhanced control and steerability. Suited for a wide range of anatomical and clinical angiography procedures.



Code	Sizes	Guide Wire Stiffness	Tip Shape
HR-2851	0.018"x150cm	Standard	Angled
HR-2051	0.018"x180cm	Standard	Angled
HR-3811	0.035"x150cm	Standard	Straight
HR-3818	0.035"x150cm	Stiff	Straight
HR-5654	0.035"x150cm	Standard	Angled
HR-5655	0.035"x150cm	Stiff	Angled
HR-6511	0.035"x180cm	Standard	Angled
HR-6512	0.035"x180cm	Standard	Straight
HR-3631	0.035"x260cm	Stiff	Angled
HR-5641	0.035"x260cm	Standard	Angled
HR-5642	0.035"x260cm	Stiff	Straight
HR-5652	0.035"x260cm	Standard	Straight
HR-4851	0.038"x150cm	Standard	Angled


Flexible Tip Length: 3 cm

Qty/Box: 5

Guide Wires

Amplatz Super Stiff

All Amplatz Guide Wires are pre-coat PTFE and can be used for interventional cardiovascular procedures and interventional biliary drainage, abscess drainage, urology procedures and catheter exchanges



Code	Sizes	Flexible Tip Length	Tip Shape
HR-3840	0.035"x80cm	10cm	Straight
HR-3815	0.035"x150cm	10cm	Straight
HR-3825	0.035"x150cm	12cm	Modified J
HR-4821	0.035"x180cm	10cm	Straight
HR-2922	0.035"x180cm	12cm	Modified J
HR-4828	0.035"x260cm	10cm	Straight
HR-3828	0.035"x260cm	12cm	Modified J
HR-4840	0.038"x80cm	10cm	Straight
HR-4815	0.038"x150cm	10cm	Straight

Qty/Box: 5

Bentson

The Bentson Guide Wire is pre-coat PTFE, with an extra floppy tip. It is used for cerebral angiography and other angiographic interventional procedures.



Code	Sizes	Flexible Tip Length	Tip Shape
HR-7912	0.035"x150cm	15cm	Straight

Qty/Box: 5

Rosen

All Rosen-Heavy Duty Guide Wires are pre-coat PTFE which offer an intermediate level of body stiffness. The non-traumatic tip characteristics are critical to the success of clinical vascular and non-vascular procedures.



Code	Sizes	Tip Shape
HR-3827	0.035"x150cm	Tight J 1.5mm
HR-3826	0.035"x180cm	Tight J 1.5mm
HR-3929	0.035"x260cm	Tight J 1.5mm

Flexible Tip Length: 1cm

Qty/Box: 5



QMD™ provides a full range of products and services to help OEMs bring medical devices to market. Specializing in silicone and other elastomerics, we develop and manufacture a broad array of high-performance single-use devices as well as precision components and subassemblies. Our end-to-end services include engineering & development, regulatory affairs, quality assurance, and supply chain support. Our global team partners with OEMs to provide tailored solutions to meet your needs.

Biometrix™

Biometrix™ offers comprehensive selection of Biometrix™ branded medical devices for endosurgery applications. Biometrix™ expertly engineers solutions and provides supply chain management for an array of clinical areas such as emergency units, operating rooms and Cath-labs.

DSP MED S.L.
QMD's Middle East Exclusive Distributor

CONTACT US >>>

DSP^{S.L.}med
Biometrix™ Products
Performance for Life

 | WWW.DSPMed.COM

 | info@dspmed.com

 | Plaza de la Solidaridad nº12, 5ª
planta 29006 Málaga
Spain