



BIOMETRIX™ PRODUCT CATALOGUE

CRITICAL CARE



Biometrix™ Critical Care Products are designed for easy access and for measuring and monitoring patients' vital signs during medical procedures performed in acute settings.

The comprehensive and diverse product portfolio has been developed to support healthcare providers with delivering efficient, patient-centered, and cost-effective care.

3 Lumen Hemodialysis Catheter Sets for Subclavian and Femoral Approach

3 Lumen Catheter Sets


Code	Description
EG-7613	12Fr x 15 cm
EG-Y7613*	12Fr x 15 cm
EG-7713	12Fr x 20 cm
EG-Y7713*	12Fr x 20 cm

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.035" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Heparin Locks (2)

(*) EG-Y7613, EG-Y7713 contains a 0.035" Nitinol Guidewire




3 Lumen Catheter Sets with Y-Introducer Needle

Code	Description
EG-I7613	12Fr x 15 cm
EG-I7713	12Fr x 20 cm

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.035" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Y Introducer Needle
- Vessel Dilators (2)
- Scalpel
- 5ml Syringe
- Heparin Locks (2)




3 Lumen Catheter Sets with Quick Insertion Syringe

Code	Description
EG-R7613	12Fr x 15 cm
EG-R7713	12Fr x 20 cm

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.035" Marked Guidewire Nitinol J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- Quick Insertion Syringe
- Heparin Locks (2)



2 Lumen Hemodialysis Catheter Sets for Subclavian and Femoral Approach


2 Lumen Catheter Sets

Code	Description
EF-7613	12Fr x 15 cm
EF-Y7613*	12Fr x 15 cm
EF-7713	12Fr x 20 cm
EF-Y7713*	12Fr x 20 cm
EF-7813	12Fr x 24 cm
EF-Y7813*	12Fr x 24 cm

Catheter Tray Contents: Qty/Box: 10

- Radiopaque PU Catheter
- 0.038" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Heparin Locks (2)

(*) EF-Y7613, EF-Y7713, EF-Y7813 contains a 0.035" Nitinol Guidewire




2 Lumen Catheter Sets with Y-Introducer Needle

Code	Description
EF-I7613	12Fr x 15 cm
EF-I7713	12Fr x 20 cm
EF-I7813	12Fr x 24 cm

Catheter Tray Contents: Qty/Box: 10

- Radiopaque PU Catheter
- 0.038" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- Scalpel
- 5ml Syringe
- Heparin Locks (2)




2 Lumen Catheter Sets with Quick Insertion Syringe

Code	Description
EF-R7613	12Fr x 15 cm
EF-R7713	12Fr x 20 cm
EF-R7813	12Fr x 24 cm

Catheter Tray Contents: Qty/Box: 10

- Radiopaque PU Catheter
- 0.035" Nitinol Guidewire J Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- Quick Insertion Syringe
- Heparin Locks (2)



3 Lumen Hemodialysis Catheter Sets for Internal Jugular Approach

3 Lumen Catheter Sets


Code	Description
EG-7623	12Fr x 15 cm curved extensions
EG-Y7623*	12Fr x 15 cm curved extensions
EG-7723	12Fr x 20 cm curved extensions
EG-Y7723*	12Fr x 20 cm curved extensions

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.035" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Heparin Locks (2)

(*) EG-Y7623, EG-Y7723 contains a 0.035" Nitinol Guidewire




3 Lumen Catheter Sets with Y-Introducer Needle

Code	Description
EG-I7623	12Fr x 15 cm curved extensions
EG-I7723	12Fr x 20 cm curved extensions

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.035" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Y Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Scalpel
- Heparin Locks (2)




3 Lumen Catheter Sets with Quick Insertion Syringe

Code	Description
EG-R7623	12Fr x 15 cm curved extensions
EG-R7723	12Fr x 20 cm curved extensions

Qty/Box: 10

Catheter Tray Contents:


- Radiopaque PU Catheter
- 0.035" Marked Nitinol Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- Quick Insertion Syringe
- Heparin Locks (2)



2 Lumen Hemodialysis Catheter Sets for Internal Jugular Approach

2 Lumen Catheter Sets

Code	Description
EF-7623	12Fr x 15 cm curved extensions
EF-Y7623*	12Fr x 15 cm curved extensions
EF-7723	12Fr x 20 cm curved extensions
EF-Y7723*	12Fr x 20 cm curved extensions



Catheter Tray Contents:


Qty/Box: 10

- Radiopaque PU Catheter
- 0.038" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Heparin Locks (2)

(*) EF-Y7623, EF-Y7723 contains a 0.035" Nitinol Guidewire

2 Lumen Catheter Sets with Y-Introducer Needle

Code	Description
EF-I7623	12Fr x 15 cm curved extensions
EF-I7723	12Fr x 20 cm curved extensions




Catheter Tray Contents:

Qty/Box: 10

- Radiopaque PU Catheter
- 0.038" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Y Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Scalpel
- Heparin Locks (2)

2 Lumen Catheter Sets with with Quick Insertion Syringe

Code	Description
EF-R7623	12Fr x 15 cm curved extensions
EF-R7723	12Fr x 20 cm curved extensions



Catheter Tray Contents:


Qty/Box: 10

- Radiopaque PU Catheter
- 0.035" Marked Nitinol Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- Quick Insertion Syringe
- Heparin Locks (2)

2 Lumen Hemodialysis Pre-Curved Catheter Sets


2 Lumen Catheter Sets	
Code	Description
EF-7693	11.5Fr x 15 cm curved extensions
EF-Y7693*	11.5Fr x 15 cm curved extensions
EF-7793	11.5Fr x 20 cm curved extensions
EF-Y7793*	11.5Fr x 20 cm curved extensions

- Catheter Tray Contents:**
- Radiopaque PU Catheter
 - 0.038" Marked Guidewire J-Flex Tip
 - 18Ga x 7cm Introducer Needle
 - Vessel Dilators (2)
 - 5ml Syringe
 - Heparin Locks (2)
- (* EF-Y7693, EF-Y7793 contains a 0.035" Nitinol Guidewire)




2 Lumen Catheter Sets with Y-Introducer Needle	
Code	Description
EF-I7693	12Fr x 15 cm curved extensions
EF-I7793	12Fr x 20 cm curved extensions

- Catheter Tray Contents:**
- Radiopaque PU Catheter
 - 0.038" Marked Guidewire J-Flex Tip
 - 18Ga x 7cm Y Introducer Needle
 - Vessel Dilators (2)
 - 5ml Syringe
 - Scalpel
 - Heparin Locks (2)



2 Lumen Catheter Sets with with Quick Insertion Syringe	
Code	Description
EF-R7693	12Fr x 15 cm curved extensions
EF-R7793	12Fr x 20 cm curved extensions

- Catheter Tray Contents:**
- Radiopaque PU Catheter
 - 0.035" Marked Guidewire Nitinol J-Flex Tip
 - 18Ga x 7cm Introducer Needle
 - Vessel Dilators (2)
 - Quick Insertion Syringe
 - Heparin Locks (2)



High Flow Hemodialysis Catheters Sets

High Flow Hemodialysis Catheter Sets specifically used for hemofiltration and rapid hemodialysis.


2 Lumen Catheter Sets	
Code	Description
EF-8613	14Fr x 15 cm
EF-F8613*	14Fr x 15 cm
EF-8713	14Fr x 20 cm
EF-F8713*	14Fr x 20 cm
EF-8813	14Fr x 24 cm
EF-F8813*	14Fr x 24 cm

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.038" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Heparin Locks (2)

(*) EF-F8613, EF-F8713, EF-F8813 contains a 0.035" Nitinol Guidewire and Adhesive Wound Dressing




2 Lumen Catheter Sets Jugular Pre-Curved	
Code	Description
EF-8693	14Fr x 15 cm
EF-F8693*	14Fr x 15 cm
EF-8793	14Fr x 20 cm
EF-F8793*	14Fr x 20 cm

Qty/Box: 10

Catheter Tray Contents:

- Radiopaque PU Catheter
- 0.038" Marked Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle
- Vessel Dilators (2)
- 5ml Syringe
- Heparin Locks (2)

(*) EF-F8693, EF-F8793 contains a 0.035" Nitinol Guidewire and Adhesive Wound Dressing





QMD™ provides a full range of products and services to help OEMs bring medical devices to market. Specializing in silicone and other elastomerics, we develop and manufacture a broad array of high-performance single-use devices as well as precision components and subassemblies. Our end-to-end services include engineering & development, regulatory affairs, quality assurance, and supply chain support. Our global team partners with OEMs to provide tailored solutions to meet your needs.

Biometrix™

Biometrix™ offers comprehensive selection of Biometrix™ branded medical devices for endosurgery applications. Biometrix™ expertly engineers solutions and provides supply chain management for an array of clinical areas such as emergency units, operating rooms and Cath-labs.

DSP MED S.L.
QMD's Middle East Exclusive Distributor

CONTACT US >>>

DSP^{S.L.}med
Biometrix™ Products
Performance for Life

 | WWW.DSPMed.COM

 | info@dspmed.com

 | Plaza de la Solidaridad nº12, 5ª
planta 29006 Málaga
Spain