



# **BIOMETRIX™ PRODUCT CATALOGUE**

# ANGIOGRAPHY & CATH-LAB

## DRAINAGE & NEPHROSTOMY CATHETERS



Interventional Cardiology and Radiology healthcare professionals rely on Biometrix™ products for a complete line of Diagnostic, Therapeutic Devices and Accessories.

Through constructive dialogues with partners and customers, we are dedicated to developing effective solutions and contributing to global medical advancement.

QMD™ offers a complete product line for interventional cardiology and radiology. Biometrix™ Contrast media injection lines, sets and syringes, angiography guide wires, haemostatic valve introducer sets, PTCA accessories kits and comprehensive custom procedure trays -part of Biometrix™ CPT program - are some of the Angio-Line™ products and accessories.

### Hemostatic Valve Introducer Sets


#### Short Length Pediatric Introducers

Code	Dilator Size	Needle	Packaging
HD-P6201	4Fr	+	Tray
HD-P6301	5Fr	+	Tray

Qty/Box: 10

##### Set Contents:

- 7cm Length Sheath/Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.018" x 40cm Double Ended Guidewire
- 21Ga x 3.8cm Introducer Needle




#### Standard Length Introducers

Code	Dilator Size	Needle	Packaging
HD-6301	5Fr	+	Tray
HD-6401	6Fr	+	Tray
HD-6501	7Fr	+	Tray
HD-6701	8Fr	+	Tray
HD-6801	9Fr	+	Tray
HD-7301	5Fr	-	Pouch
HD-7401	6Fr	-	Pouch
HD-7501	7Fr	-	Pouch
HD-7701	8Fr	-	Pouch
HD-7901	9Fr	-	Pouch

Qty/Box: 10

##### Set Contents:

- 12cm Length Sheath/Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.038" x 50cm Double Ended Guidewire
- 18Ga x 7cm Introducer Needle




#### Long Length Introducers

Code	Dilator Size	Packaging
HD-9301	5Fr	Pouch
HD-9401	6Fr	Pouch
HD-9501	7Fr	Pouch
HD-9701	8Fr	Pouch

Qty/Box: 10

##### Set Contents:

- 23cm Length Sheath/Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.038" x 70cm Double Ended Guidewire
- 18Ga x 7cm Introducer Needle




Micropuncture Hemostatic Valve Introducer Sets

Standard Length Introducers		
Code	Dilator Size	Packaging
HD-Y6201	4Fr	Tray
HD-Y6301	5Fr	Tray
HD-Y6401	6Fr	Tray

Qty/Box: 10

Set Contents:

- 12cm Length Sheath/Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 21Ga x 3.8cm Introducer Needle
- 0.018" x 40cm Stainless Steel Guidewire
- 5ml Syringe




Long Length Introducers		
Code	Dilator Size	Packaging
HD-X9301	5Fr	Pouch
HD-X9401	6Fr	Pouch
HD-X9501	7Fr	Pouch

Qty/Box: 10

Set Contents:

- 23cm Length Sheath/Dilator Hemostatic Valve Unit with Side Arm and Stopcock
- 0.018" x 70cm Stainless Steel Guidewire
- 21Ga x 3.8cm Introducer Needle



Split Sheath Introducer Sets	
Code	Description
JF-0701	7Fr
JF-0801	8Fr
JF-0901	9Fr
JF-1001	10Fr
JF-1101	11Fr

Qty/Box: 10

Set Contents:


- Dilator Assembly
- 18Ga x 7cm Introducer Needle
- 0.038" x 50cm Double Ended Guidewire
- 5ml Syringe L/L



Biometrix™ Micropuncture Sets offer a smooth and atraumatic entry into the patient's vasculature. The 0.018" Nitinol Guidewire with its platinum tip allows vascular access which is achieved in a few short steps with minimal vessel trauma.

### Micropuncture Sets

Micropuncture Sets		
Code	Description	Guide Wire Material
HT-0210	4Fr	Nitinol
HT-0310	5Fr	Nitinol



Qty/Box: 10


#### Set Contents:

- 4Fr or 5Fr Vessel Dilator Set
- 0.018" x 40cm Nitinol Guidewire
- 21 Ga x 7cm Introducer Needle with Ecogenic Tip

### Advantages of Micropuncture


- Minimizes trauma and potential complications (hematoma)
- Minimizes venous spasm which can complicate or hinder the procedure
- Reduces post procedure mechanical phlebitis

Guidewires					
Guidewire Code	Diameter	Length	Tip	J Radius	Material
HR-2242	0.018"	40cm	Platinum Floppy 3cm	3mm	Nitinol
HR-2542	0.018"	65cm	Platinum Floppy 3cm	3mm	Nitinol



Qty/Box: 50

Introducer Needles				
Code	Cannula Size	No. of parts	Bevel	Wing
HQ-1510	18Ga x 7cm	One Part	Standard	-
HQ-1511	18Ga x 7cm	One Part	Standard	+
HQ-1514	18Ga x 7cm	One Part	Short	-
HQ-1521	18Ga x 7cm	Two Part	Standard	-
HQ-3510	18Ga x 7cm	Y-Needle	Standard	-
HQ-1525	18Ga x 9cm	One Part	Standard	-
HQ-1611	19Ga x 7cm	One Part	Standard	+
HQ-1710	20Ga x 4cm	One Part	Standard	+
HQ-1810	21Ga x 3.8cm	One Part	Standard	-



Qty/Box: 100



**QMD™** provides a full range of products and services to help OEMs bring medical devices to market. Specializing in silicone and other elastomerics, we develop and manufacture a broad array of high-performance single-use devices as well as precision components and subassemblies. Our end-to-end services include engineering & development, regulatory affairs, quality assurance, and supply chain support. Our global team partners with OEMs to provide tailored solutions to meet your needs.

### **Biometrix™**

Biometrix™ offers comprehensive selection of Biometrix™ branded medical devices for endosurgery applications. Biometrix™ expertly engineers solutions and provides supply chain management for an array of clinical areas such as emergency units, operating rooms and Cath-labs.


**DSP MED S.L.**  
QMD's Middle East Exclusive Distributor

**CONTACT US >>>**

**DSP<sup>S.L.</sup>med**  
**Biometrix™ Products**  
Performance for Life

 | [WWW.DSPMed.COM](http://WWW.DSPMed.COM)

 | [info@dspmed.com](mailto:info@dspmed.com)

 | Plaza de la Solidaridad nº12, 5ª  
planta 29006 Málaga  
Spain