



BIOMETRIX™ PRODUCT CATALOGUE

ANGIOGRAPHY & CATH-LAB

DRAINAGE & NEPHROSTOMY CATHETERS



Interventional Cardiology and Radiology healthcare professionals rely on Biometrix™ products for a complete line of Diagnostic, Therapeutic Devices and Accessories.

Through constructive dialogues with partners and customers, we are dedicated to developing effective solutions and contributing to global medical advancement.

Manifold & Contrast Media Sets & Angiography Waste Bags

Contrast Media Manifold Sets		
Code	Description	Qty/Box
BQ-1000	2 Port Angiography Manifold Set	25

Set Contents:

- Manifold 2 port F/F with 15cm High Pressure Line with rotating adaptor
- Infusion Set Vented
- 12ml Finger Grip Syringe

Code	Description	Qty/Box
BQ-2000	Contrast Media Set	25

Set Contents:

- 150cm Infusion Line
- Double Check Valve
- 3 Way High Pressure Stopcock
- Infusion Set Vented
- 12ml Finger Grip Syringe

Code	Description	Qty/Box
BQ-3000	3 Port Angiography Manifold Set	25






Set Contents:

- 3 Port Manifold F/MLL Rot.
- 100cm Infusion Line F/MLL
- Infusion Set Vented
- 12ml Finger Grip Syringe

Angiography Waste Bag

The Angiography waste bag provides safe and efficient contrast media management.

Code	Description	Qty/Box
BQ-0205	Double-Line Angiography Waste Bag	25
BQ-0215	Single-Line Angiography Waste Bag	25





QMD™ provides a full range of products and services to help OEMs bring medical devices to market. Specializing in silicone and other elastomerics, we develop and manufacture a broad array of high-performance single-use devices as well as precision components and subassemblies. Our end-to-end services include engineering & development, regulatory affairs, quality assurance, and supply chain support. Our global team partners with OEMs to provide tailored solutions to meet your needs.

Biometrix™

Biometrix™ offers comprehensive selection of Biometrix™ branded medical devices for endosurgery applications. Biometrix™ expertly engineers solutions and provides supply chain management for an array of clinical areas such as emergency units, operating rooms and Cath-labs.


DSP MED S.L.
QMD's Middle East Exclusive Distributor

CONTACT US >>>

DSP^{S.L.}med
Biometrix™ Products
Performance for Life

 | WWW.DSPMed.COM

 | info@dspmed.com

 | Plaza de la Solidaridad nº12, 5ª
planta 29006 Málaga
Spain