



BIOMETRIX™ PRODUCT CATALOGUE

SURGERY



Biometrix™ is a leading supplier of surgical essential devices. From chest and wound drainage systems to lumbar puncture sets, the products are produced to high manufacturing standards and exact quality specifications.

PNEUMOTHORAX SETS


Biometrix™ Pneumothorax sets are used for removing both air and fluid from the pleural space. Available in a variety of sizes and design configurations, provide quick set up and safe treatment solutions.

Code	Description
NP-1000	Non Seldinger Type Pneumothorax Set

Qty/Box: 10

Set Contents:

- Radiopaque Drainage Catheter with 6Fr Dilator Assembly Needle (20Ga)
- Connecting Tube with 3 Way Stopcock
- Scalpel
- 10ml Syringe
- Stepped Adaptor




Code	Description
NP-1001	Seldinger Type Pneumothorax Set

Qty/Box: 10

Set Contents:

- Radiopaque Drainage Catheter with 6Fr Dilator Assembly
- Connecting Tube with 3 Way Stopcock
- 0.032" Stiff Guidewire J-Flex Tip
- 10ml Syringe
- 18Ga x 7cm Introducer Needle
- Scalpel #11 with Short Handle
- Stepped Adaptor
- 800ml Drainage Bag
- Drape 50x50cm
- Heimlich Valve




Code	Description
NP-1005	Seldinger Type Pneumothorax Set

Qty/Box: 10

Set Contents:

- Radiopaque Drainage Catheter with 6Fr Dilator Assembly
- Connecting Tube with 3 Way Stopcock
- 0.032" Stiff Guidewire J-Flex Tip
- 18Ga x 7cm Introducer Needle Echogenic Tip
- Scalpel
- 10ml Syringe
- Stepped Adaptor
- Dilator 11Fr x 15cm




PIGTAIL PNEUMOTHORAX SETS

Pneumothorax Pigtail Catheter Set One Step Technique


Code	Description
NP-2000	14Fr Pneumothorax Pigtail Catheter Set

Set Contents:

- 14Fr X 30cm Pneumothorax Catheter with straightener
- Metal Stiffening Cannula
- Metal Trocar Needle
- Connecting Tube with 3 way stopcock for sampling option & closing
- Scalpel
- Heimlich Valve



Qty/Box: 5




Pneumothorax Pigtail Catheter Set Seldinger Technique


Code	Description
NP-2001	14Fr Pneumothorax Pigtail Catheter Set Seldinger Technique

Set Contents:

- 14Fr X 30cm Pneumothorax Catheter with straightener
- Metal Stiffening Cannula
- Flexible Stiffening Cannula
- Connecting Tube with 3 way stopcock for sampling option & closing
- 0.038" X 80cm Amplatz Guide wire
- 18Ga X 9cm Introducer Needle
- Vessel Dilators
- Scalpel
- Heimlich Valve



Qty/Box: 5




Pediatric Pneumothorax Pigtail Catheter Set Seldinger Technique


Code	Description
NP-2100	8Fr Pediatric Pneumothorax Pigtail Catheter Set Seldinger Technique

Set Contents:

- 18Fr X 20cm Pneumothorax Catheter with safety string lock
- Connecting Tube with 3 way stopcock for sampling option and closing
- 0.028 X 40cm J Amplatz Guide wire
- 18Ga X 7cm Introducer Needle
- 8Fr Vessel Dilator
- Scalpel
- Heimlich Valve



Qty/Box: 10





QMD™ provides a full range of products and services to help OEMs bring medical devices to market. Specializing in silicone and other elastomerics, we develop and manufacture a broad array of high-performance single-use devices as well as precision components and subassemblies. Our end-to-end services include engineering & development, regulatory affairs, quality assurance, and supply chain support. Our global team partners with OEMs to provide tailored solutions to meet your needs.

Biometrix™

Biometrix™ offers comprehensive selection of Biometrix™ branded medical devices for endosurgery applications. Biometrix™ expertly engineers solutions and provides supply chain management for an array of clinical areas such as emergency units, operating rooms and Cath-labs.


DSP MED S.L.
QMD's Middle East Exclusive Distributor

CONTACT US >>>

DSP^{S.L.}med
Biometrix™ Products
Performance for Life

 | WWW.DSPMed.COM

 | info@dspmed.com

 | Plaza de la Solidaridad nº12, 5ª
planta 29006 Málaga
Spain